

# Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US2004/014325

International filing date: 07 May 2004 (07.05.2004)

Document type: Certified copy of priority document

Document details: Country/Office: US  
Number: 60/469,383  
Filing date: 09 May 2003 (09.05.2003)

Date of receipt at the International Bureau: 16 January 2007 (16.01.2007)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland  
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

1562748

# THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

*January 05, 2007*

**THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.**

**APPLICATION NUMBER: 60/469,383**

**FILING DATE: *May 09, 2003***

**RELATED PCT APPLICATION NUMBER: *PCT/US04/14325***

**THE COUNTRY CODE AND NUMBER OF YOUR PRIORITY APPLICATION, TO BE USED FOR FILING ABROAD UNDER THE PARIS CONVENTION, IS *US60/469,383***



Certified by

Under Secretary of Commerce  
for Intellectual Property  
and Director of the United States  
Patent and Trademark Office

05/09/03  
JCS39 U.S. PTO

05-12-03

60469383 .05/09/03

PTO/SB/16 (10-01)

Approved for use through 10/31/2002. OMB 0651-0032

U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

# **PROVISIONAL APPLICATION FOR PATENT COVER SHEET**

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No. EY 013434076 US

PTO  
J1040 U.S. PTO  
60/469383

05/09/03

INVENTOR(S)					
Given Name (first and middle (if any))		Family Name or Surname		Residence (City and either State or Foreign Country)	
Robert C. Jaime		Getts Schwalm		Collegeville, Pennsylvania Norristown, Pennsylvania	
<input type="checkbox"/> Additional inventors are being named on the _____ separately numbered sheets attached hereto					
TITLE OF THE INVENTION (500 characters max)					
Methods for Amplification of Nucleic Acid Sequences					
Direct all correspondence to:			CORRESPONDENCE ADDRESS		
<input checked="" type="checkbox"/> Customer Number <u>28083</u>			<div style="border: 1px solid black; padding: 5px; text-align: center;">             Place Customer Number Bar Code Here           </div>		
OR			Type Customer Number here		
<input type="checkbox"/> Firm or Individual Name			<b>Morris E. Cohen</b> <small>PATENT TRADEMARK OFFICE</small>		
Address			1122 Coney Island Avenue		
Address			Suite 217		
City		State	NY	ZIP	11230
Country		US	Telephone	718-859-8009	Fax 718-859-3044
ENCLOSED APPLICATION PARTS (check all that apply)					
<input checked="" type="checkbox"/> Specification		Number of Pages <u>33</u>		<input type="checkbox"/> CD(s), Number <u>        </u>	
<input checked="" type="checkbox"/> Drawing(s)		Number of Sheets <u>10</u>		<input type="checkbox"/> Other (specify) <u>        </u>	
<input type="checkbox"/> Application Data Sheet. See 37 CFR 1.76					
METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT					
<input type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.			FILING FEE AMOUNT (\$)		
<input checked="" type="checkbox"/> A check or money order is enclosed to cover the filing fees			<div style="border: 1px solid black; padding: 10px; text-align: center;"> <b>\$160.00</b> </div>		
<input checked="" type="checkbox"/> The Commissioner is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number <u>50-1604</u>					
<input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.					
The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.					
<input type="checkbox"/> No.					
<input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____					

Respectfully submitted,

SIGNATURE *M. Cohen*

TYPED or PRINTED NAME Morris E. Cohen

TELEPHONE 718-859-8009

Date 05/09/2003

REGISTRATION NO. 39,947

(if appropriate)

Docket Number: 4081.025

## **USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT**

This collection of information is required by 37 CFR 1.51. The information is used by the public to file (and by the PTO to process) a provisional application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the complete provisional application to the PTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, Washington, D.C. 20231. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Box Provisional Application, Assistant Commissioner for Patents, Washington, D.C.

## Methods for Amplification of Nucleic Acid Sequences

**Inventors:** Robert C. Getts      Collegeville, Pennsylvania  
                  Jaime Schwalm      Norristown, Pennsylvania

### Background of the Invention

Methods for the amplification of nucleic acid sequences are currently important tools in molecular biology. Frequently, a nucleic acid is present in extremely minute quantities in a given sample, making it advantageous to provide techniques for amplification of the nucleic acid to sufficient amounts to facilitate its further analysis.

Several such amplification techniques are currently employed in the art. In one well known method, for example, disclosed by Eberwine and his colleagues, a compound oligonucleotide is utilized for the amplification, wherein the compound oligonucleotide is provided with both a T7 promoter and a primer. (See e.g. U.S. Patent Nos. 5,716,785; 5,891,636; 5,958,688; and 6,291,170 B1; all of which are fully incorporated herein by reference).

Briefly summarized, in the Eberwine method, a cDNA copy is created of an initial mRNA, with subsequent second strand synthesis to create a cDNA which is double stranded. Amplification is conducted via the promoter from the compound oligonucleotide, with transcription proceeding off of the cDNA's second strand. Since the second strand is used for transcription, the Eberwine method produces amplified RNA which is anti-sense to the initial mRNA sequence. As a result, it produces materials which are the wrong sense for use on oligonucleotide arrays, when using conventional labeling



methods (i.e. reverse transcription and dye incorporation) and not incorporating the label into the aRNA.

Furthermore, since transcription occurs in Eberwine in the 3' to 5' direction with respect to the initial mRNA sequence, the method is, by its nature, biased toward that 3' end of the initial mRNA sequence, the end where priming occurs and from where transcription begins. In the Eberwine method, if the second strand template is not completely intact, a bias can be introduced into the results, favoring messages at the 3' end of the initial mRNA sequence, with loss of messages at the 5' end. For example, if second strand synthesis does not proceed to completion, messages corresponding to the 5' end of the initial mRNA will not be transcribed, resulting in an incomplete copy of the initial message. Thus, in general, the method may result in amplification wherein not every message is duplicated during the amplification procedure.

Accordingly, it would be advantageous to provide an improved method in the art for amplification of nucleic acid sequences. It would be particularly advantageous to provide a method which produces amplified nucleic acid of the same sense as the initial nucleic acid sample. It would further be advantageous to provide a method of amplification in which the copying of the amplified sequences results in a messages having complete fidelity to both the 5' and 3' ends of the initial nucleic acid sequence.



Further in accordance with the invention, methods are provided for amplification of nucleic acid wherein priming occurs from one end of the nucleic acid strand and transcription occurs from the other end of the strand, to ensure the highest fidelity possible between the initial nucleic acid strand and the amplified copy. In the preferred embodiments, priming is conducted from the 3' end of the initial nucleic acid with a promoter subsequently added to a sequence corresponding to the 5' end of the initial nucleic acid, such that priming and nucleic acid synthesis initially occurs in the 3' to 5' direction, subsequently followed by transcription in the 5' to 3' direction relative to the initial nucleic acid sequence of interest.

In accordance with one embodiment of the preferred method, an existing mRNA sequence is used to create a single stranded cDNA. For example, in one preferred embodiment, a dT primer can be hybridized to the mRNA's 3' poly A tail and extended to form a double stranded mRNA/cDNA molecule, followed by degradation of the mRNA strand. Alternatively, any other suitable method for creation of a single stranded cDNA from the mRNA can be used. For example a random primer or a gene specific primer can be used to form the double stranded mRNA/cDNA molecule, followed by degradation of the mRNA.

Upon creation of the single stranded cDNA molecule (ss cDNA), a dA tail is attached to the molecule using any suitable protocol. In the preferred embodiment, a dA tail is attached to the 3' end using the enzyme Terminal deoxynucleotidyl Transferase (TdT).

Alternatively, a tail of another specifically desired sequence can be added to the 3' end. In the event that a tail other than a dA tail is added, the subsequent steps using a dT sequence complementary to the dA tail (for hybridization to that tail), are modified to include the appropriate sequence instead of dT, i.e. to use the sequence complementary to the tail that has been attached. Or, in the event that the sequence at the 5' end of the initial nucleic acid is known, addition of the tail to the 3' end can be skipped, if desired. When the sequence at the 5' end is known, the complement to that sequence can be treated as equivalent to the tail for further steps; the complement being present at the 3' end of the cDNA in the case of RNA as the initial nucleic acid, for example. However, addition of a tail (whether dA or otherwise) is preferred and will generally be necessary due to the fact that, in most applications, the sequence of the 5' end of the initial nucleic acid is unknown.

As an alternative to the use of mRNA, a strand of DNA can be used as the initial nucleic acid from the sample. In this embodiment, a random primer is used to create a complementary strand to the initial strand of DNA (the complementary strand being analogous to the cDNA described above), with a tail then being added to the complementary strand (e.g. through attachment to the 3' end by TdT) for use in the same manner as the dA tailed cDNA described herein (or with use of an existing sequence as the tail, as previously discussed).

Once the tailed nucleic acid has been produced (or a molecule is used with an existing sequence used as the tail), a promoter is placed on the 3' end by any of various suitable methods. The promoter is used to initiate transcription off of the 3' end of the tailed molecule. In other words, transcription proceeds off of the dA sequence, in the example when a dA tail is used. This transcription, proceeding

in the 3' to 5' direction with respect to the molecule with the promoter, corresponds to copying in the 5' to 3' direction with respect to the initial nucleic acid of interest. By the use of multiple rounds of transcription from the promoter, large numbers of copies are created of the initial nucleic acid, the copies being referred to as amplified mRNA or "Amp-mRNA" herein (when the initial nucleic acid is RNA), or as amplified DNA (when the initial nucleic acid is DNA).

Each of the amplified copies is identical in sense to the initial nucleic acid message. For example, in the case of an initial mRNA, since the multiple transcripts/amplified copies are all complementary to the cDNA, which is itself complementary to the initial mRNA, the amplified mRNA copies all contain the identical message to the initial mRNA sequence; i.e. amplified sense copies of the initial mRNA are produced, and not anti-sense copies.

These sense copies of amplified nucleic acid can then be used in any manner in which nucleic acids are used in the art, including any current, past, or future protocols.

In the case of amplified mRNA, in preferred embodiments, the amplified copies are converted to amplified cDNA copies, which can then be used in desired assays. In further preferred embodiments, the amplified nucleic acids are used in assays using microarrays. In alternative or additional embodiments of the invention, the amplified nucleic acids are used in labeled with dendritic reagents, particularly dendrimers such as 3DNA<sup>TM</sup> reagents available from Genisphere, Inc. of Montvale, New Jersey. Dendritic reagents are further described in Nilsen et al., Dendritic Nucleic Acid Structures, J. Theor. Biol., 187, 273-284 (1997); in Stears et al., A Novel, Sensitive Detection System for High-

Density Microarrays Using Dendrimer Technology, *Physiol. Genomics*, 3: 93-99 (2000); and in various U.S. patents, such as U.S. Patent Nos. 5,175,270; 5,484,904; 5,487,973; 6,072,043; 6,110,687; and 6,117,631, these references all being fully incorporated herein by reference.

### **Brief Description of the Figures**

Figure 1 is a schematic illustration of a preferred method for amplification of nucleic acids in accordance with present invention, using staggered ligation to a tailed cDNA followed by nucleic acid amplification.

Figure 2 is a schematic illustration of an additional preferred method for amplification of nucleic acids in accordance with the present invention, using blunt ligation to a tailed cDNA, followed by nucleic acid amplification.

Figure 3 is a schematic illustration of a further additional preferred method for amplification of nucleic acids using hybridization of a complex primer to a tailed cDNA followed by nucleic acid amplification.

### **Detailed Description of the Invention and the Preferred Embodiments**

In accordance with the present invention, methods are provided for the amplification of nucleic acid sequences. In the preferred embodiments, methods are provided for the amplification of DNA or RNA to produce amplified sense copies of the initial sequence. Further in the preferred embodiments, methods are provided wherein the amplification step (i.e. the step in which the amplified copies are

generated) utilizes nucleic acid copying commencing from the end corresponding to the 5' end of the initial sequence, via transcription in the 5' to 3' direction with respect to the initial nucleic acid strand.

In the preferred embodiments, a dA tailed cDNA copy is produced from the initial RNA strand. (Alternatively, a strand of DNA can serve as the initial nucleic acid, with a complement being created that is tailed and then used in the same manner as the dA tailed cDNA described herein). A promoter is then placed on the 3' dA tail, with the promoter being used to initiate multiple rounds of transcription to produce amplified sense copies of the initial nucleic acid, with transcription proceeding off of the dA tailed strand. Thus, as opposed to the Eberwine method, whether or not the second strand is intact will not impact the fidelity of the amplified transcript, so that loss of message is minimized. In fact, second strand synthesis is optional in the procedures of the present invention, as shown for example, in Figure 1 attached hereto, so that the present invention is effective even if the second strand is not synthesized. In addition, the methods of the present invention require significantly less time than the Eberwine method. Whereas, the Eberwine protocol normally takes several days, the present process can be done in one day if desired, or two if it is desired to conduct it overnight.

The amplified DNA or RNA can then be used in any desired context or protocol. Where the amplified nucleic acid is amplified mRNA, direct use of the amplified mRNA in further assays can be conducted for example, following such methods as disclosed in PCT Application Serial No. PCT/US01/22818 filed July 19, 2001 (Int'l Publication No. WO 02/06511), which is fully incorporated herein by reference.

004-695077-0050003

The foregoing discussion discloses and describes merely exemplary embodiments of the present invention. One skilled in the art will readily recognize from such discussion, and from the accompanying drawings and examples, that various modifications and variations can be made therein without departing from the spirit and scope of the invention. The present application is intended to cover all such modifications and variations.





having the initial mRNA hybridized to a complementary cDNA. (Or, as discussed above, a random primer or gene specific primer can be used in place of the dT primer, to create the molecule of mRNA/cDNA). The mRNA strand of this duplex is then degraded, leaving a single stranded cDNA molecule ("nucleic acid c"). In the case of use of a dT primer, the single stranded cDNA ("ss cDNA") has a 5' end having a poly dT or dU tail, the 5' dT (or dU) tail having originated from the primer. Likewise, in the case of the use of another primer, that primer will remain at the 5' end of the single stranded cDNA.

In step 3, this single stranded cDNA is modified to add a dA tail thereto (this step also being referred to herein as "tailing"). In the preferred embodiment, the enzyme Terminal deoxynucleotidyl transferase (TdT) is used to add a dA tail to the 3' terminus of the cDNA ("TdT/dA tailing"), TdT being a known DNA polymerase catalyzing the terminal addition of deoxyribonucleotides to the 3'-OH termini of DNA. At the end of this step, a cDNA molecule has been produced ("nucleic acid d") which is complementary to the initial mRNA sequence, and which has a 3' dA tail attached thereto.

In step 4 of the method shown in Figure 1, a promoter is attached to the dA tailed cDNA's 3' end. In a preferred embodiment, the T7 promoter is utilized, although other suitable promoters can alternatively be used consistent with the present invention. In the method of Figure 1, staggered ligation is conducted of a double stranded T7 promoter to the 3' end of the tailed cDNA. Alternatively, other methods may be used for attachment of the promoter to the tailed cDNA, including, but not limited to those described below, and shown in Figures 2 and 3.

In the method of staggered ligation shown in Figure 1 (step 4), a double stranded (“ds”) promoter is provided having a dT overhang, with the dT overhang hybridizing the double stranded T7 promoter to the cDNA’s 3’ dA tail. The 3’ dA tail of the dA tailed cDNA molecule, therefore, provides a partial overlapping complementary sequence to the ds T7 promoter, thereby aligning the T7 promoter for ligation. Using a ligase, the double stranded T7 promoter is ligated to the dA tailed cDNA molecule, so that the dA tailed cDNA now carries a promoter on its 3’ end (“nucleic acid e”).

Following ligation of the promoter onto the dA tailed cDNA, a large number of mRNA messages can be transcribed off of the cDNA, as shown in Figure 1, step 5. For example, when the T7 promoter is used, T7 polymerase and NTPs can be added to transcribe messages off of the cDNA template, as shown in step 5 of Figure 1. In the method of Figure 1, transcription can be conducted off of a single stranded cDNA provided that the double stranded promoter is affixed thereto. Alternately, in alternative embodiments, transcription can be conducted from a double stranded molecule, as shown, for example, in the method of Figure 2.

Thus, in step 5 (the amplification step), repeated cycles of transcription are conducted to create large numbers of copies of the initial mRNA in the sample, the copies being referred to herein as amplified mRNA (“Amp-mRNA”; “nucleic acid f” in Figure 1). The amplified mRNA copies include a sequence identical to the mRNA initially present in the sample and, furthermore, are sense copies of the initial mRNA, not antisense. Furthermore, transcription is conducted in the 3’ to 5’ direction with respect to the dA tailed cDNA, which corresponds to copying in the 5’ to 3’ direction with respect to the initial mRNA sequence.





These amplified cDNAs can then be used in a desired experiment or assay. Any desired assay means can be used. In the preferred embodiment, the assay uses a microarray for simultaneously conducting a large number of large scale "test-tube" like experiments.

Likewise, any desired labeling means can be used for detection of the cDNAs in the assay. For example, the cDNAs can be produced using radioactive nucleotides, for labeling and use of the cDNAs in subsequent experiments.

In a preferred embodiment, as shown in step 8 of Figure 1, a capture sequence is ligated onto the end of each amplified cDNA, the capture sequence being complementary to a label molecule. Examples of such labeling methods using capture sequences are disclosed in the PCT Applications previously cited herein and cited below, all of which are fully incorporated herein by reference.

When the cDNA is applied to the assay medium (e.g. to a microarray surface), the label molecule having the complement to the capture sequence is used to detect the hybridization of the cDNA to the array. Any desired molecule can be used for the label. In preferred embodiments, the label molecule is a labeled dendrimer to provide increased labeling efficiency. Dendrimers and the use of dendrimers is described in Nilsen et al., Dendritic Nucleic Acid Structures, *J. Theor. Biol.*, 187, 273-284 (1997); in Stears et al., A Novel, Sensitive Detection System for High-Density Microarrays Using Dendrimer Technology, *Physiol. Genomics*, 3: 93-99 (2000); and in various patents and patent applications, such as U.S. Patent Nos. 5,175,270; 5,484,904; 5,487,973; 6,072,043; 6,110,687; 6,117,631, PCT Application Serial Number PCT/US01/07477 filed March 8, 2001 (Int'l Publication

No. WO 01/066555), PCT Application Serial No. PCT/US01/22818 filed July 19, 2001 (Int'l Publication No. WO 02/06511), and PCT Application Serial No. PCT/US01/29589 filed September 20, 2001 (Int'l Publication No. WO 02/033125), all of which are fully incorporated herein by reference.

In a further preferred embodiment of the invention, Locked Nucleic Acid ("LNA") is used as a blocker to prevent non-specific hybridizations, as further discussed in PCT Application Serial No. PCT/US02/27799 filed 3 September 2002 (Int'l Publication No. WO 03/020902), which is fully incorporated herein by reference. For example, as shown in step 9 of Figure 1, the LNA blocker can be used to hybridize to the addition located at the 3' end of the amplified cDNA to minimize any nonspecific binding effects when the amplified cDNA is applied to the assay medium.

Any desired method can then be used to conduct the desired assay. For example, in one preferred embodiment, a two step hybridization can be conducted as shown in step 10 of Figure 1. In a two step hybridization, hybridization of the cDNA to the array and hybridization of the label molecule to the cDNA (via the capture sequence) are conducted as two separate steps. Alternatively, a one step hybridization can be conducted as discussed in PCT Application Serial Number PCT/US01/07477 filed 8 March 2001 (Int'l Publication No. WO 01/066555), which is fully incorporated herein by reference.

A hybridization pattern is then generated, and signal detection is conducted to view and then analyze the assay results.

## *Method 2 - Blunt Ligation of T7 Promoter*

In accordance with an alternative method for amplification of nucleic acid sequences, blunt ligation of a T7 promoter to a dA tailed cDNA is conducted, as shown for example, in Figure 2.

In this method, a user begins with a nucleic acid sequence (preferably mRNA), and using that sequence generates a 3' dA tailed cDNA, in the same manner as previously described above with respect to method 1.

Following creation of the dA tailed cDNA, a dT primer is hybridized onto the 3' dA tail, as shown in step 4 of Figure 2. Using DNA polymerase and dNTPs, a double stranded cDNA is then generated ("nucleic acid f" of Figure 2), the double stranded molecule being provided with a blunt end at the end corresponding to the 3' dA tailed cDNA, as a result of the exonuclease activity of the DNA polymerase.

Using DNA ligase, a double stranded T7 promoter is then ligated onto the double stranded cDNA (producing "nucleic acid g"). The promoter is ligated onto the end of the ds cDNA to which the 3' dA tail had been added, as shown in Step 5 of Figure 2.

T7 polymerase and NTPs are then used to transcribe mRNA messages from the double stranded cDNA, producing amplified mRNA ("nucleic acid h" of Figure 2). In other words, amplified sense copies of the original mRNA are produced, with transcription having been conducted in a direction corresponding to the 5' to 3' direction of the initial mRNA.



The amplified mRNAs, can now be used in any desired method and for any desired purpose, as previously discussed. In preferred embodiments, amplified cDNAs can be created, e.g. for subsequent microarray assays as discussed with respect to Figure 1; or, alternatively, any other desired protocol can be conducted.

### *Method 3 - Complex Primer - Alternative T7 Amplification*

In accordance with a further alternative method for amplification of nucleic acid sequences, a method for amplification using a complex primer is conducted, as shown for example, in Figure 3.

In this method, a user begins with a nucleic acid sequence (preferably mRNA), and using that sequence generates a 3' dA tailed cDNA, in generally the same manner as previously described above. As shown in Figure 3, a dT primer can be used as previously discussed. Alternatively, in this step a random primer or gene specific primer can be used, as discussed in conjunction with Method 1. Likewise, a random primer or gene specific primer can be alternately utilized in method 2 described above, as well.

Once the dA tailed cDNA has been produced, a single strand of a T7 promoter having a dT tail is hybridized to the 3' dA end of the cDNA, as shown in step 4 of Figure 3. Using DNA polymerase and dNTPs, the single stranded cDNA with attached single stranded promoter is converted to a double stranded molecule, including conversion of the single stranded T7 promoter to a double stranded T7 promoter.

Subsequently, in step 5, using T7 polymerase and NTPs, amplified mRNA is produced from the double stranded cDNA molecule. The amplified mRNA is of the same sense as the initial mRNA, and is produced using transcription in a direction corresponding to the 5' to 3' direction of the initial mRNA strand. As discussed above, this amplified mRNA can be used in any method or context in which mRNA is utilized in molecular biological protocols. For example, the amplified mRNA can be used with the methods shown in steps 6 onward of Figure 3, and/or in any of the methods described or illustrated with respect to Methods 1 and 2 (and corresponding Figures 1 and 2).

Having described the invention with respect to several preferred embodiments, several illustrative examples of protocols using the present inventions are provided below.

#### Example 1:

##### Opposite Strand T7-Amplification (Oligo Array Compatible) using Staggered Ligation

###### First Strand Synthesis:

- a. Reverse transcribe up to 10ug of total RNA using cap03 Genisphere RT primer in a total volume of 20ul. (Superscript II or Genisphere RT Reagents may be used). Note: For RNA amounts greater than 5ug, use 5pmol/ul RT primer. Standard dT and/or RP primers and/or standard RT primers may be used instead of Genisphere primers.

1-10ul total RNA

1ul 1pmol/ul or 5pmole/ul Genisphere RT primer

---

11ul

Heat to 80°C for 10 minutes. Ice immediately for 1-2 minutes.

Add RT components using one of the following options:

4ul SSII 5X 1 <sup>st</sup> Strand Buffer	4ul 5X Rxn Buffer (Promega)
2ul DTT	2ul nuclease free water (Vial 10)
1ul dNTP (Genisphere 10mM ea)	1ul dNTP mix (Vial 3)
1ul Supersasin	1ul RNase Inhibitor (Vial 4)
1ul Superscript II	1ul MMLV RT (Promega)
<hr/>	
20ul final volume	

Incubate at 42°C for 1.5-2 hours.

- b. Stop reaction with 3.5ul of 0.5M NaOH/50mM EDTA and heat to 65°C for 10 minutes.
- c. Neutralize the reaction with 5ul of 1M Tris-HCl, pH 7.5.
- d. Adjust volume to 50ul with 1X TE pH 8.0.
- e. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.

#### TdT Tailing of 1<sup>st</sup> Strand cDNA:

##### Protocol Including Qiagen Purification of TdT Reaction:

- a. Heat cDNA to 80°C for 5-10 minutes to denature. Cool on ice immediately for 1-2 minutes.

- b. Tail cDNA with dATP mix to about 20-50bp length as described below, and incubate the reaction at 37°C for 10-60 minutes (TBD).

50ul	cDNA
14ul	5X TdT buffer
2ul	dATP mix, no ddATP (~10mM dATP stock)
2.5ul	water
1.5ul	TdT
<hr/>	
70ul	

**Alternatively:** Tris based buffer may be used in place of standard TdT buffer  
avoid purification prior to ligation (see Protocol Excluding Purification, on following page).

- c. Stop the reaction with 8ul 0.5M EDTA. Add 22ul 1X TE pH 8.0 to adjust the volume to 100ul.
- d. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.
- e. Heat to 95-100°C for 10-15 minutes to denature cDNA, ice immediately for 1-2 minutes.
- f. Add T7 promoter ligation components as follows and incubate reaction at room temperature for 30-60 minutes.

50ul	tailed cDNA
10ul	6X Ligation Mix (containing T7 promoter sequence)
2.5ul	T4 DNA Ligase
<hr/>	
62.5ul	

- g. Add 7ul 0.5M EDTA. Optional: Heat to 80°C for 10minutes. Add 30.5ul of 1X TE pH 8.0 to adjust final volume to 100ul.
- h. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.

**Protocol Excluding Qiagen Purification of TdT Reaction:**

- a. Heat cDNA to 80°C for 5-10 minutes to denature. Cool on ice immediately for 1-2 minutes.
- b. Tail cDNA with dATP mix to about 20-50bp length as described below, and incubate the reaction at 37°C for 10-60 minutes (TBD).

50ul cDNA
6ul 10X Tris-based buffer
2ul dATP mix, no ddATP (~10mM dATP stock)
1ul water
1ul TdT
<hr/>
60ul

- c. Incubate at 95-100°C for 10-15 to inactivate TdT enzyme and denature cDNA, ice immediately for 1-2 minutes. (DO NOT USE EDTA!)
- d. Add ligation components as follows and incubate reaction at room temperature for 30-60 minutes.

60.0ul tailed cDNA
12.5ul 6X Ligation Mix
2.5ul T4 DNA Ligase
<hr/>
75.0ul

- e. Add 8ul 0.5M EDTA. Optional: Heat to 80°C for 10minutes. Add 17ul of 1X TE pH 8.0 to adjust final volume to 100ul.
- f. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.

**Concentration of cDNA Containing T7 Promoter Sequences:**

- g. Ethanol precipitate cDNA to concentrate.
- h. Resuspend pellet in nuclease free water. Heat to 80°C for 15-20 minutes, vortex.
- i. Remove up to 2ug equivalent for IVT reaction. Note: When using Ambion messageAmp aRNA kit, aliquot must not exceed 8ul.

**In Vitro Transcription:**

- a. Add in vitro transcription reagents and incubate reaction at 37°C for 6-14 hours (see below).

8ul purified cDNA (containing ligated T7 promoter)  
 2ul T7 ATP Soln (75mM)  
 2ul T7 CTP Soln (75mM)  
 2ul T7 GTP Soln (75mM)  
 2ul T7 UTP Soln (75mM)  
 2ul T7 10X Reaction Buffer  
 2ul T7 Enzyme Mix (including RNase inhibitor)

---

20ul (use Ambion reagents or equivalent)

- b. Optional: Add 2ul DNase I to reaction and incubate for 30 minutes at 37°C.
- c. Optional: Heat to 80°C for 10 minutes to inactivate DNase I. Note: DNase I will be removed during purification.

#### **Purification of Amplified RNA Using Qiagen RNeasy Kit**

- a. Purify amplified RNA sample using Qiagen RNeasy kit (protocol for RNA cleanup) as directed by the manufacturer.
- b. QC RNA sample by obtaining O.D. reading and visualizing sample on gel.

#### **Example 2: Opposite Strand T7-Amplification (Oligo Array Compatible) Using Blunt End Ligation**

##### **First Strand Synthesis:**

- a. Reverse transcribe up to 10ug of total RNA using cap03 Genisphere RT primer in a total volume of 20ul. (Superscript II or Ambion messageAMP aRNA Reagents may be used).
- Note: For RNA amounts greater than 5ug, use 5pmol/ul RT primer. Standard dT and/or RP primers and/or standard RT primers may be used instead of Genisphere primers.

1-10ul total RNA	
1ul 1pmol/ul or 5pmole/ul Genisphere RT primer	
11ul	

Heat to 80°C for 10 minutes. Ice immediately for 1-2 minutes.

Add RT components using one of the following options:

4ul SSII 5X 1 <sup>st</sup> Strand Buffer	2ul 10X 1 <sup>st</sup> Strand Buffer (Ambion)
2ul DTT	1ul nuclease free water (Ambion)
1ul dNTP (Genisphere 10mM ea)	1ul dNTP mix (Ambion)
1ul Supersasin	1ul RNase Inhibitor (Ambion)
1ul Superscript II	1ul MMLV RT (Ambion)
<hr/>	
20ul final volume	

Incubate at 42°C for 1.5-2 hours.

- b. Stop reaction with 3.5ul of 0.5M NaOH/50mM EDTA and heat to 65°C for 10 minutes.
- c. Neutralize the reaction with 5ul of 1M Tris-HCl, pH 7.5.
- d. Adjust volume to 50ul with 1X TE pH 8.0.
- e. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.

#### TdT Tailing of 1<sup>st</sup> Strand cDNA:

- a. Heat cDNA to 80°C for 5-10 minutes to denature. Cool on ice immediately for 1-2 minutes.
- b. Tail cDNA with dATP mix to about 50bp length as described below, and incubate the reaction at 37°C for 30-60 minutes.



50ul cDNA  
 14ul 5X TdT buffer  
 2ul dATP mix (50:1 ratio of 100mM dATP:ddATP)  
 2.5ul water  
 1.5ul TdT

---

70ul

- c. Stop the reaction with 8ul 0.5M EDTA. Add 22ul 1X TE pH 8.0 to adjust the volume to 100ul.
- d. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.

#### Second Strand Synthesis:

- a. Add 2ul standard dT primer (about 100 pmol) to the purified first strand reaction. Heat the mixture to 95-100°C for 10 minutes to denature possible A/T circular hybrids. Ice immediately for 1-2 minutes.
- b. Add remaining second strand reagents (see below) and incubate the reaction at 16°C for 2 hours.

50ul purified cDNA  
 22ul water  
 10ul 10X second strand buffer  
 2ul dT primer (added above)  
 4ul dNTP mix  
 2ul DNA polymerase I

---

100ul (use Ambion reagents or other equivalent)



[illegible]

- 8ul purified 2<sup>nd</sup> strand cDNA  
2ul T7 ATP Soln (75mM)  
2ul T7 CTP Soln (75mM)  
2ul T7 GTP Soln (75mM)  
2ul T7 UTP Soln (75mM)  
2ul T7 10X Reaction Buffer  
2ul T7 Enzyme Mix

- b. Optional: Add 2ul DNase I to reaction and incubate for 30 minutes at 37°C.
- c. Optional: Heat to 80°C for 10 minutes to inactivate DNase I. Note: DNase I will be removed during purification.

- a. Purify amplified RNA sample using Qiagen RNeasy kit (protocol for RNA cleanup) as directed by the manufacturer.
- b. QC RNA sample by obtaining O.D. reading and visualizing sample on gel.

**Example 3:  
Opposite Strand T7-Amplification (Oligo Array Compatible)  
Using T7(dT) 24 Primer**

**First Strand Synthesis:**

- a. Reverse transcribe up to 10ug of total RNA using cap03 Genisphere RT primer in a total volume of 20ul. (Superscript II or Ambion messageAMP aRNA Reagents may be used).

Note: For RNA amounts greater than 5ug, use 5pmol/ul RT primer. Standard dT and/or RP primers and/or standard RT primers may be used instead of Genisphere primers.

1-10ul total RNA
1ul 1pmol/ul or 5pmole/ul Genisphere RT primer
<hr/>
11ul

Heat to 80°C for 10 minutes. Ice immediately for 1-2minutes.

Add RT components using one of the following options:

4ul SSII 5X 1 <sup>st</sup> Strand Buffer	2ul 10X 1 <sup>st</sup> Strand Buffer (Ambion)
2ul DTT	1ul nuclease free water (Ambion)
1ul dNTP (Genisphere 10mM ea)	1ul dNTP mix (Ambion)
1ul Supersasin	1ul RNase Inhibitor (Ambion)
1ul Superscript II	1ul MMLV RT (Ambion)
<hr/>	
20ul final volume	

Incubate at 42°C for 1.5-2 hours.

- b. Stop reaction with 3.5ul of 0.5M NaOH/50mM EDTA and heat to 65°C for 10 minutes.
- c. Neutralize the reaction with 5ul of 1M Tris-HCl, pH 7.5.
- d. Adjust volume to 50ul with 1X TE pH 8.0.
- e. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.

#### TdT Tailing of 1<sup>st</sup> Strand cDNA:

- a. Heat cDNA to 80°C for 5-10 minutes to denature. Cool on ice immediately for 1-2 minutes.
- b. Tail cDNA with dATP mix to about 50bp length as described below, and incubate the reaction at 37°C for 30-60 minutes.

50ul cDNA
14ul 5X TdT buffer
2ul dATP mix (50:1 ratio of 100mM dATP:ddATP)
2.5ul water
1.5ul TdT
<hr/>
70ul

- c. Stop the reaction with 8ul 0.5M EDTA. Add 22ul 1X TE pH 8.0 to adjust the volume to 100ul.
- d. Purify the reaction using Qiagen QIAquick PCR purification kit as directed by the manufacturers protocol. Elution will be in 50ul.

## Second Strand Synthesis:

- a. Add 2ul of T7-dT(24) primer (about 100pmol) to the purified first strand reaction.

Heat the mixture to 95-100°C for 10 minutes to denature possible A/T circular hybrids. Ice immediately for 1-2 minutes.

- b. Add remaining second strand reagents (see below) and incubate the reaction at 16°C for 2 hours.

50ul purified cDNA
22ul water
10ul 10X second strand buffer
2ul T7-dT(24) primer (added above)
4ul dNTP mix
2ul DNA polymerase I
<hr/>
100ul (use Ambion reagents or other equivalent)

- c. Purify second strand reaction using Qiagen QIAquick PCR Purification kit as directed by the manufacturers protocol. Elution will be in 50ul.
- d. Ethanol precipitate cDNA to concentrate.
- e. Resuspend pellet in nuclease free water. Heat to 80°C for 15-20 minutes, vortex.
- f. Remove up to 2ug equivalent for IVT reaction. Note: When using Ambion messageAmp aRNA kit, aliquot must not exceed 8ul.

### In Vitro Transcription:

- a. Add in vitro transcription reagents and incubate reaction at 37°C for 6-14 hours  
(see below).

8ul purified 2<sup>nd</sup> strand cDNA  
 2ul T7 ATP Soln (75mM)  
 2ul T7 CTP Soln (75mM)  
 2ul T7 GTP Soln (75mM)  
 2ul T7 UTP Soln (75mM)  
 2ul T7 10X Reaction Buffer  
 2ul T7 Enzyme Mix

---

20ul (use Ambion reagents or equivalent)

- b. Optional: Add 2ul DNase I to reaction and incubate for 30 minutes at 37°C.
- c. Optional: Heat to 80°C for 10 minutes to inactivate DNase I. Note: DNase I will be removed during purification.

### Purification of Amplified RNA Using Qiagen RNeasy Kit

- a. Purify amplified RNA sample using Qiagen RNeasy kit (protocol for RNA cleanup) as directed by the manufacturer.
- b. QC RNA sample by obtaining O.D. reading and visualizing sample on gel.

The foregoing discussion, therefore, discloses and describes merely exemplary and preferred embodiments of the present invention. One skilled in the art will readily recognize from such discussion, and from the accompanying drawings and examples, that various changes, modifications and variations can be made therein without departing from the spirit and scope of the invention as

defined in the claims. One skilled in the art may likewise by applying current or future knowledge, adopt the same for use in accordance with the present invention. Yet, having described this invention with regard to specific embodiments, it is to be understood that the description is not meant as a limitation, and that the present application covers all such embodiments, modifications and variations of the presently disclosed inventions.


### Claims

What is claimed is:

1. Methods for amplification of nucleic acid sequences having all of the features disclosed above.
2. Nucleic acid sequences produced by any of the methods for amplification of nucleic acid sequences disclosed above.



Figure 1: Staggered Ligation to tailed cDNA  
followed by T7 Amplification



1 dA tailed cDNA (nucleic acid D)



dsT7 promoter

Ligase (4)

dT sequence



(5)

T7 Polymerase  
& NTPs

sequence  
from promoter  
dT sequence



(Nucleic Acid F)

Reverse  
Transcriptase

(6)

Array 350  
Random Primer

(Nucleic Acid G)

Amp-cDNA

Amp-mRNA




Degrade RNA (7)

A-cDNA  
(Nucleic Acid H)




(PAGE 35)





## Figure 2: Blunt Ligation of T7 Promotor

for use with Genisphere Array 350 or Array350RP Kits  
Initial dT based priming follow by random prime target labeling



dA tailed cDNA 3' 5' } NUCLEIC ACID E

dT primer  
DNA Polymerase & dNTPs

(4)

3' 5' } NUCLEIC ACID F  
3' 5' }

(5)

DNA Ligase

ds T7 promoter

ds T7 promoter


3' 5' } NUCLEIC ACID G  
3' 5' }

(6)


T7 Polymerase & NTPs

Amp-mRNA 5' 3'  
(NUCLEIC ACID H)

(PAGE 38)







# Figure 3: Complex Primer Alternative T7 Amplification

for use with Genisphere Array 350 ~~(dTT primer)~~ <sup>(dT primer) kit</sup> Poly A or dA

